

Protein Kinase Filter Binding Assay Protocol (Manual Assay)

Standard Procedure for all kinases

Materials:

ATP	(e.g. Sigma, # A-7699)
CaCl₂	(e.g. VWR, # 1.02382)
Calmodulin	(e.g. Millipore, # 14-368)
cGMP	(e.g. Sigma, # G-6129)
1,2-Dioleoyl-glycerol	(e.g. Sigma, # D-8394)
DNA	(e.g. Sigma, # D-4522)
DMSO	(e.g. Sigma, # 154938)
DTT	(e.g. Sigma, # D-0632)
EDTA	(e.g. Roth, # 8040.1)
H₂O p.a.	(e.g. VWR, # 1.16754.9010)
HEPES	(e.g. Sigma, # H-3375)
MgCl₂	(e.g. Sigma, # M-3634)
MnCl₂	(e.g. VWR, # 1.05927.1000)
NaCl	(e.g. Merck, # 1.06404)
Na-orthovanadate	(e.g. Sigma, # S-6508)
Ortho-phosphoric acid	(e.g. VWR, # 1.00563.1000)
PEG_{20.000}	(e.g. SERVA, # 33138)
Phosphatidyl-serine	(e.g. Fluka, # 79405)
(10 mg/ml in chloroform/methanol)	
³³P-γ-ATP	(e.g. PerkinElmer, # NEG 302 H)
96-well V-bottom MTP	
(Polypropylen)	
	(e.g. Nunc, # 442587)
96-well Multiscreen Vacuum Manifold	
(e.g. Millipore, # MAVM096OR)	
96-well Multiscreen Filterplates	
(depending on used substrate):	
- Phosphocellulose Membrane	
	(e.g. Millipore, # MSPHN6B10)
- GlassFiber Membrane	
	(e.g. Millipore, # MSFCN6B10)
Scintillation Cocktail	(e.g. Fluka, # 03999)

Prepare in advance (store at 4° C):

Standard-Assay-Buffer (2.5x stock)

- 125 mM HEPES-NaOH, pH 7.5
- 7.5 mM MgCl₂
- 7.5 mM MnCl₂
- 7.5 μM Na-orthovanadate
- 2.5 mM DTT

Kinase-Dilution-Buffer (10x stock)

- 500 mM HEPES-NaOH, pH7.5
- 2.5 mg/ml PEG_{20.000}
- 10 mM DTT

Final Assay Conc. (50 μl reaction cocktail):

- 70 mM HEPES-NaOH, pH 7.5
- 3 mM MgCl₂
- 3 mM MnCl₂
- 3 μM Na-orthovanadate
- 1.2 mM DTT
- 50 μg/ml PEG_{20.000}
- ATP [var. conc. (dep. on K_M) + ³³P-γ-ATP]
- 1 % (v/v) DMSO
- Substrate (variable)
- Recombinant Protein Kinase (variable)

All PKC assays (except the PKC-mu and the PKC-nu assay) additionally contain 1 mM CaCl₂, 4 mM EDTA, 5 μg/ml Phosphatidylserine and 1 μg/ml 1,2-Dioleoyl-glycerol.

The CAMK1D, CAMK2A, CAMK2B, CAMK2D, CAMK4, CAMKK1, CAMKK2, DAPK2, EEF2K, MYLK, MYLK2 and MYLK3 assays additionally contain 1 μg/ml Calmodulin and 0.5 mM CaCl₂.

The PRKG1 and PRKG2 assays additionally contain 1 μM cGMP.

The DNA-PK assay additionally contains 2.5 μg/ml DNA.

Manual Assay Procedure:

In PP V-bottom plate:

1. Add 20 μl Standard-Assay-Buffer (2.5x stock)
2. Add 5 μl 10% DMSO (*w, w/o test cpd*)
3. Add 10 μl substrate (in 50mM HEPES pH 7.5)
4. Add 10 μl recombinant protein kinase (in 1x Kinase-Dilution-Buffer)
5. Add 5 μl ATP (in H₂O)
6. Mix on shaker
7. Incubate for 40 min at 30°C
8. Stop reaction with 20 μl 10% H₃PO₄
9. Prewet filterplate with 200 μl 150 mM H₃PO₄
(wash MSPH-filterplate with 20 μl 100% EtOH first)
10. **Transfer reaction mix to filterplate**
11. Incubate for at least 10, but no longer than 20 min
12. Filtrate reaction mix by applying vacuum
13. Wash 3 times with 200 μl 150 mM H₃PO₄
14. Wash once with 20 μl 100% EtOH
15. Dry for 30 min at 40°C
16. Add 50 μl scintillator per well
17. Count with Scintillation Counter