## HybriDetect

Universal lateral flow dipstick for the detection of biotin- and FITC-labeled analytes (proteins, genomic amplificates);  
development platform  
English: Page 1-9

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**Note:** Significant changes are indicated by dotted lines in the margin. In the end of the IFU you will find a table listing the reasons for the introduced changes.

## Table of Contents

Explanation of Symbols .....	3
Warnings and Precautions .....	3
Intended Use .....	3
Materials Supplied, Storage and Stability .....	4
Materials Required (not included) .....	4
Necessary Development Work - Development Platform .....	4
Method .....	5
Common Testprinciple .....	6
Control Band Dipstick .....	6
Interpretation of Results .....	7
Assay Performance "PCR-Products" .....	8
Interpretation of „PCR“ Results .....	8
Trouble Shooting „PCR“ .....	9
Assay Sensitivity „PCR“ .....	9
Literature References .....	9
Additional Available Products .....	9

## Explanation of Symbols

Symbols (GB) Symbole (DE)	Explanation Erklärung	Symbols (GB) Symbole (DE)	Explanation Erklärung
	Expiry date Haltbarkeitsdatum		Package size Packungsgröße
	Consult Instructions for Use Gebrauchsanweisung beachten		Manufacturer Hersteller
	Batch code Los-Bezeichnung		Storage conditions Lagerungsbedingungen
	Catalogue number Katalognummer		Consult attended documents Begleitdokumente beachten

## Warnings and Precautions

- All reagents should be stored at 2 – 8° C in their original containers.
- Before use, bring all reagents to room temperature (18 – 28° C).
- The expiration date of all components must be observed.
- Protect dipsticks from humidity; Container must always be closed.
- Touch and label only the foil- covered areas of the dipsticks (labeling area).
- The disposal of waste materials must be carried out according to current local regulations.
- The assay buffer contains an anti-microbial reagent; therefore avoid contact with skin and/or mucous membranes.
- For professional users

## Intended Use

HybriDetect is a universal lateral flow dipstick for the detection of analytes (proteins, genomic amplicates) labelled with FITC and biotin. It is a development platform. The test is for research use only, not for diagnostic purposes.

## Materials Supplied, Storage and Stability

Components	Cat.-No.	Content	Preparation	Storage	Shelf Life
<b>HybriDetect Dipsticks:</b> Membrane coated with biotin-ligand; polyclonal (rabbit) anti-FITC antibody in gold conjugate	MGDS	2 x 50 tests	Ready-to use	2 –8°C; Container must be always closed (protected from moisture)!	until expiry date
<b>HybriDetect Assay Buffer:</b> Tris-buffered saline	MGCB	2 vials à 10 mL	Ready-to use	2 –8°C	until expiry date

Material Safety Data Sheets are available on request or from the company website ([www.milenia-biotec.de](http://www.milenia-biotec.de))

## Materials Required (not included)

- Pipets
- Pipet tips (containing protective filters for PCR)
- Reaction tubes or 96-well microtiter plate

## Necessary Development Work - Development Platform

Development of a solution, containing two different labeled detectors for the analyte.

Following conditions are necessary:

- 1) Detectors must be labeled with:
  - FITC (fluorescein isothiocyanate)
  - Biotin
- 2) Use about 100 µL fluid (sample material and analyte-specific solution) for the assay procedure.
- 3) The provided assay buffer (MGCB) may be used as a basic buffer for this analyte-specific solution.

Example of use:

20 µl sample material and 80 µl analyte-specific solution with 5 min incubation time.

**Notice:**

Volumes, analyte-specific solution, and incubation time are part of the individual test development.

A basic procedure for the detection of genomic amplicons is explained on page 7: "Assay performance PCR products". The amplification product should be biotinylated and the specific hybridization probe should be FITC labeled. In principle all amplification procedures (Polymerase Chain Reaction or isothermal amplifications like LAMP or RPA) can be used.

## Method


Milenia<sup>®</sup> HybriDetect is a ready-to-use, universal test strip (dipstick), which is based on the lateral flow technology using gold particles. The dipstick is designed to develop qualitative or semi-quantitative rapid test systems for several analytes such as proteins, antibodies, or gene amplificons. The user needs to develop an analyte-specific solution, which contains a first detector (e.g. antibody, antigen, specific probe) labeled with FITC and a second one (e.g. antibody, primer) labeled with biotin, (see "Common Testprinciple" on page 5).

The sample to be determined is mixed with the developed analyte-specific solution, and then the dipstick is placed into the solution.

The complexed analyte, labeled with FITC and biotin, binds first to the gold-labeled FITC-specific antibodies in the sample application area of the dipstick. The gold complexes travel through the membrane, driven by capillary forces. Only the analyte captured gold particles will bind when they pass the line with the immobilized biotin-ligand molecules and generate a red-blue band over the time. Unbound gold particles migrate over the control band and will be captured by species-specific antibodies. With prolonged incubation time, the formation of an intensely colored control band appears.

## Common Testprinciple


1) Ready-to-use universal dipstick (analyte-independent)

### Symbols:


■	FITC
●	Biotin
●	Goldparticle
●	anti-FITC Antibody
●	Analyte Detector 1 + 2
●	Analyte (Protein, antibody, genetic amplicon)
■	Biotin-Ligand
Y	anti-rabbit Antibody

2) Addition of sample and analyte-specific solution

3) Incubation

4) Result:

## Control Band Dipstick


In any case, the control band must be visible!

It is a control function and can not be used to assess the quality of the result of the test band. If the control band is not visible after the incubation period, the result is invalid! The test must be repeated with a new dipstick!


## Interpretation of Results

There are two possibilities to interpret Milenia<sup>®</sup> HybriDetect dipstick results:

### 1) Qualitative by visual interpretation


### 2) Semi-quantitative, e. g. by evaluation cards

## Assay Performance “PCR-Products”

1. Take the required number of dipsticks out of the container and mark them.
2. For each sample to be analyzed pipet 100 µL HybriDetect Assay Buffer or individual developed buffer into a reaction tube or a well of a microtiter plate.
3. Pipet 5 – 10 µl of the hybridization product directly on the sample application area, or alternatively: add 5 – 10 µl of the hybridization product into the solution of the reaction tube / well.
4. Place the dipsticks with the sample application area into the solution and incubate them e. g. for 5 - 15 minutes in an upright position.
5. In the end of incubation period, remove dipsticks from assay solution and interpret test results immediately.

### Notice:

- If a higher analytical sensitivity is required, it could be helpful to increase the volume of the PCR product.
- Volumes, analyte-specific solution, and incubation time are always part of the individual test development.

## Interpretation of „PCR“ Results

<b>1. Test band <u>and</u> control band are clearly formed red bands</b> <b>Notes:</b> 1.1 Weakly stained test bands have to be regarded as positive 1.2 Positive results may be visible within 1-2 minutes. 1.3 In case of very high concentrations of hybridization product, control band's intensity may go down. Nevertheless, the control band should be still visible clearly.	The PCR amplicon was detected <b>(positive)</b> .
<b>2. Only control band is visible as red line.</b>	<b>No</b> PCR amplicon was detected <b>(negative)</b> .

positive      negative      not valid

A diagram showing three dipstick results. The first dipstick on the left has two distinct red bands and is labeled 'positive'. The second dipstick in the middle has one distinct red band and is labeled 'negative'. The third dipstick on the right has no visible red bands and is labeled 'not valid'. To the right of the dipsticks, two arrows point to the bands. The top arrow is labeled 'Control Band' and the bottom arrow is labeled 'Test Band'.

### Important note:

To check specificity of the reaction the introduction of a PCR negative control is strongly recommended.

## Trouble Shooting „PCR“

Problem	Possible cause(s)	Recommendation
Control band is not visible.	a) wrong or degraded assay buffer b) expiration date of dipsticks is exceeded c) wrong storage conditions of dipsticks	apply new (fresh) chemicals
Negative result with dipstick but clearly visible band in agarose gel	a) detection of an unspecific PCR product in agarose gel b) hybridization was not successful	Check identity of PCR product by Southern blotting oder sequence analysis Check conditions of hybridization.
Mineral oil	a) mineral oil affects flow characteristics of the assay b) development of test band might be hampered.	Remove PCR product very slowly from the bottom of the reaction vial.

## Assay Sensitivity „PCR“

Analytical sensitivity of the Milenia® HybriDetect is equivalent or up to 1.000 fold higher than a Ethidium Bromide stained agarose gel. Dependent on the size of the PCR product and the number of amplification cycles up to 5 pg DNA can be detected.

## Literature References

- Kiatpathomchai W, et al, J Virol Methods (2008); 153: 214-217  
Puthawibool T, et al, J Virol Methods (2009); 156 (1-2): 27-31  
Jaroenram W, et al, Mol Cell Probes (2009); 23 (2): 65-70  
Nimitphak T, et al, Mol Cell Probes (2009); 24 (1): 1-5  
Kikuchi T, et al, Nematology (2009) ; 99 (12): 1365-1369  
Piepenburg O, et al, PLoS Biology 2006, Volume 4, Issue 7, e204

## Additional Available Products

Product Name	Order No.	Content	Description
Milenia® HybriDetect 2T	MGHD2 1	100 tests	Dipsticks with two different test bands (as third labeling with digoxigenine is required)

For any questions and suggestions, please contact us:

**Telefon:** +49 641 – 94 88 83 – 0  
**Fax:** +49 641 – 94 88 83 – 80  
**Web:** [www.milenia-biotec.de](http://www.milenia-biotec.de)

## HybriDetect

Universal lateral flow dipstick for detection of biotin- and FITC-labeled analytes (proteins, genomic amplificates);  
development platform  
English: Page: 1-9

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