


***Data*sheet**

Biotin-PEG3-Azide



Product number	Amount
BCFA-021-1	1 mg
BCFA-021-5	5 mg

Molecular formula: C₁₈H₃₂N₆O₅S

Molecular weight: 444.55 g/mol

Storage conditions: -20 °C, dry, inert gas

Purity: ≥ 90% (HPLC)

Appearance: off-white to pale brown amorphous solid

Solubility: DMSO, DMF, MeOH

Shelf life: 12 months

CAS: 875770-34-6

For research use only!