

Certificate of Analysis

Teknova Controlled Document:CofA

| | |
|-----------|----------|
| DVC-Rev: | 10674-02 |
| Rev By: | JLA |
| Check By: | BT |

Product Name / Specification: 20X PBS with 1% Tween-20. 1000mL, Sterile.

Version:

5.2

Part Number: P1192


Date QC Released: 2/26/2020

LotNumber: P119224B2001

Released By:

Date Manufactured: 2/21/2020

Date of Re-assay: 2/21/2023



Bella Neufeld, QC Manager

Appearance: Clear, colorless liquid

Method of Sterilization: 0.2µm Sterile filtration

Storage Conditions: Room temperature, 2-30°C.

QA Testing:

Result:

Expected:

Conductivity

181.0mS/cm

N/A

pH

6.33

N/A

pH1x

7.22

7.20 ± 0.04

Sterility

Pass

Pass