



### Uniform Polystyrene Latex Beads

All Uniform Polystyrene Latex Beads are in aqueous suspension at 10% solids. The average diameter tolerance will be within 10% of each catalog number specified. Product comes in 5 mL, 15 mL, and 100 mL quantity packages.. For bulk quantity of 500+ mL price, or diameter not listed here, please contact Magsphere Inc.

Catalog No.	Diameter (um)	5 mL	15 mL	100 mL	500 + mL
PS020NM	0.02	\$75	\$135	\$500	Bulk Rate
PS030NM	0.03	\$75	\$135	\$500	Bulk Rate
PS040NM	0.04	\$75	\$135	\$500	Bulk Rate
PS050NM	0.05	\$75	\$135	\$500	Bulk Rate
PS060NM	0.06	\$75	\$135	\$500	Bulk Rate
PS080NM	0.08	\$75	\$135	\$500	Bulk Rate
PS100NM	0.1	\$75	\$135	\$500	Bulk Rate
PS150NM	0.15	\$75	\$135	\$500	Bulk Rate
PS200NM	0.2	\$75	\$135	\$500	Bulk Rate
PS300NM	0.3	\$75	\$135	\$500	Bulk Rate
PS500NM	0.5	\$75	\$135	\$500	Bulk Rate
PS700NM	0.7	\$75	\$135	\$500	Bulk Rate
PS900NM	0.9	\$75	\$135	\$500	Bulk Rate
PS002UM	2	\$95	\$180	\$810	Bulk Rate
PS003UM	3	\$105	\$200	\$900	Bulk Rate
PS004UM	4	\$105	\$200	\$900	Bulk Rate
PS005UM	5	\$105	\$200	\$900	Bulk Rate
PS007UM	7	\$125	\$250	\$1,125	Bulk Rate
PS010UM	10	\$155	\$300	\$1,350	Bulk Rate

### NIST Traceable Size Standards

Our NIST traceable size standards, Monosized Polystyrene or PMMA Microspheres, are in aqueous suspension at 2% solids. Those NIST traceable beads have a CV (Coefficient of Variation) 1% to 8% for larger than 2 um, and 5% to 15% CV for smaller than 2 um particles and submicron particles. Product will be shipped within 5% of the diameter you selected. We offer in 10 mL, 20 mL, and 100 mL quantity packages. For bulk quantity of 500+ mL please contact Magsphere Inc.

Catalog No.	Diameter (um)	10 mL	20 mL	100 mL	500 + mL
NST100NM	0.1	\$125	\$180	\$600	Bulk Rate
NST200NM	0.2	\$125	\$180	\$600	Bulk Rate
NST300NM	0.3	\$125	\$180	\$600	Bulk Rate
NST500NM	0.5	\$125	\$180	\$600	Bulk Rate
NST750NM	0.75	\$125	\$180	\$600	Bulk Rate
NST001UM	1	\$155	\$225	\$750	Bulk Rate
NST002UM	2	\$155	\$225	\$750	Bulk Rate
NST003UM	3	\$155	\$225	\$750	Bulk Rate
NST004UM	4	\$155	\$225	\$750	Bulk Rate
NST005UM	5	\$155	\$225	\$750	Bulk Rate
NST006UM	6	\$155	\$225	\$750	Bulk Rate
NST007UM	7	\$155	\$225	\$750	Bulk Rate
NST008UM	8	\$155	\$225	\$750	Bulk Rate
NST010UM	10	\$155	\$225	\$750	Bulk Rate
NST015UM	15	\$155	\$225	\$750	Bulk Rate
NST020UM	20	\$155	\$225	\$750	Bulk Rate
NST030UM	30	\$155	\$225	\$750	Bulk Rate
NST040UM	40	\$155	\$225	\$750	Bulk Rate

### How to Order:

Telephone: 1-626-584-9446 Fax:1-626-795-0797

www.magsphere.com

1993 E. Locust Street, Pasadena, CA 91107

### Carboxylated Polystyrene Latex Beads

All Uniform Carboxylated Polystyrene Latex Beads are in aqueous suspension at 10% solids. The average diameter tolerance will be within 10% of each catalog number specified. Product comes in 5 mL, 15 mL, and 100 mL quantity packages.. For bulk quantity of 500+ mL price, or diameter not listed here, please contact Magsphere Inc.

Catalog No.	Diameter (um)	5 mL	15 mL	100 mL	500 + mL
CA020NM	0.02	\$85	\$160	\$580	Bulk Rate
CA040NM	0.04	\$85	\$160	\$580	Bulk Rate
CA060NM	0.06	\$85	\$160	\$580	Bulk Rate
CA080NM	0.08	\$85	\$160	\$580	Bulk Rate
CA100NM	0.1	\$85	\$160	\$580	Bulk Rate
CA200NM	0.2	\$85	\$160	\$580	Bulk Rate
CA300NM	0.3	\$85	\$160	\$580	Bulk Rate
CA400NM	0.4	\$85	\$160	\$580	Bulk Rate
CA500NM	0.5	\$85	\$160	\$580	Bulk Rate
CA750NM	0.75	\$85	\$160	\$580	Bulk Rate
CA001UM	1	\$104	\$190	\$800	Bulk Rate
CA002UM	2	\$133	\$250	\$1,120	Bulk Rate
CA005UM	5	\$185	\$350	\$1,460	Bulk Rate
CA007UM	7	\$190	\$370	\$1,750	Bulk Rate
CA010UM	10	\$190	\$370	\$1,750	Bulk Rate

### Aminated Latex Beads

All Uniform Aminated Polystyrene Latex Beads are in aqueous suspension at 10% solids. The average diameter tolerance will be within 10% of each catalog number specified. Product comes in 5 mL, 15 mL, and 100 mL quantity packages.. For bulk quantity of 500+ mL price, or diameter not listed here, please contact Magsphere Inc.

Catalog No.	Diameter (um)	5 mL	15 mL	100 mL	500 + mL
AM050NM	0.05	\$100	\$180	\$600	Bulk Rate
AM060NM	0.06	\$100	\$180	\$600	Bulk Rate
AM075NM	0.075	\$100	\$180	\$600	Bulk Rate
AM090NM	0.9	\$100	\$180	\$600	Bulk Rate
AM100NM	0.1	\$100	\$180	\$600	Bulk Rate
AM200NM	0.2	\$100	\$180	\$600	Bulk Rate
AM300NM	0.3	\$100	\$180	\$600	Bulk Rate
AM400NM	0.4	\$100	\$180	\$600	Bulk Rate
AM500NM	0.5	\$100	\$180	\$600	Bulk Rate
AM600NM	0.6	\$100	\$180	\$600	Bulk Rate
AM700NM	0.7	\$100	\$180	\$600	Bulk Rate
AM800NM	0.8	\$100	\$180	\$600	Bulk Rate
AM900NM	0.9	\$100	\$180	\$600	Bulk Rate
AM001UM	1.5	\$130	\$230	\$1,150	Bulk Rate
AM002UM	2	\$190	\$370	\$1,750	Bulk Rate

*Quality manufacturing since 1985*

Effective Date: January 10, 2012. Prices are subject to change without notice.

DVB Crosslinked Polystyrene Latex Nano-Particles and Microspheres

Polystyrene crosslinked with divinyl benzene (DVB) nano or microbeads are ideal for applications where solvent resistance or compression is needed. This product line is also in aqueous suspension at 10% solids, and the average diameter will be within 10% of the diameter catalog number specifies.

Catalog No.	Diameter (um)	5 mL	15 mL	100 mL	500 + mL
PSDVB040NM	0.04	\$85	\$160	\$580	Bulk Rate
PSDVB050NM	0.05	\$85	\$160	\$580	Bulk Rate
PSDVB060NM	0.06	\$85	\$160	\$580	Bulk Rate
PSDVB070NM	0.07	\$85	\$160	\$580	Bulk Rate
PSDVB100NM	0.1	\$85	\$160	\$580	Bulk Rate
PSDVB200NM	0.2	\$85	\$160	\$580	Bulk Rate
PSDVB400NM	0.4	\$85	\$160	\$580	Bulk Rate
PSDVB800NM	0.8	\$85	\$160	\$580	Bulk Rate
PSDVB001UM	1	\$85	\$160	\$600	Bulk Rate
PSDVB002UM	2	\$85	\$160	\$600	Bulk Rate
PSDVB003UM	3	\$90	\$170	\$650	Bulk Rate
PSDVB004UM	4	\$95	\$185	\$750	Bulk Rate
PSDVB005UM	5	\$100	\$200	\$850	Bulk Rate
PSDVB006UM	6	\$110	\$250	\$950	Bulk Rate
PSDVB007UM	7	\$125	\$300	\$1,150	Bulk Rate

PMMA Latex Beads

PMMA (Polymethyl Methacrylate) nano or microbeads are in aqueous suspension at 10% solids.Monosized polymethyl methacrylate PMMA microbeads are extra smooth and with clean surface. Product will be shipped within 10% of the diameter you selected.

Catalog No.	Diameter (um)	5 mL	15 mL	100 mL	500 + mL
PM060NM	0.06	\$80	\$150	\$525	Bulk Rate
PM070NM	0.07	\$80	\$150	\$525	Bulk Rate
PM100NM	0.1	\$80	\$150	\$525	Bulk Rate
PM150NM	0.15	\$80	\$150	\$525	Bulk Rate
PM200NM	0.2	\$80	\$150	\$525	Bulk Rate
PM300NM	0.3	\$80	\$150	\$525	Bulk Rate
PM400NM	0.4	\$80	\$150	\$525	Bulk Rate
PM600NM	0.6	\$80	\$150	\$525	Bulk Rate
PM700NM	0.7	\$80	\$150	\$525	Bulk Rate
PM800NM	0.8	\$80	\$150	\$525	Bulk Rate
PM001UM	1	\$85	\$160	\$600	Bulk Rate
PM002UM	2	\$85	\$160	\$600	Bulk Rate
PM003UM	3	\$85	\$160	\$600	Bulk Rate
PM004UM	4	\$85	\$160	\$600	Bulk Rate
PM005UM	5	\$85	\$160	\$600	Bulk Rate
PM010UM	10	\$85	\$160	\$600	Bulk Rate
PM020UM	20	\$85	\$160	\$600	Bulk Rate
PM030UM	30	\$85	\$160	\$600	Bulk Rate
PM040UM	40	\$85	\$160	\$600	Bulk Rate

Color Particles

Blue Polystyrene Latex Beads

We also offer colored polystyrene latex beads, which is also in aqueous suspension at 10% solids. You can also contact us if you are looking for a specific color intensity. Product will be shipped within 10% of the diameter you specify. Our PSB300NM and PSB400NM are widely used in rapid tests.

Catalog No.	Diameter (um)	5 mL	15 mL	100 mL	500 + mL
PSB030NM	0.03	\$103	\$195	\$850	Bulk Rate
PSB050NM	0.05	\$103	\$195	\$850	Bulk Rate
PSB060NM	0.06	\$103	\$195	\$850	Bulk Rate
PSB080NM	0.08	\$103	\$195	\$850	Bulk Rate
PSB100NM	0.1	\$103	\$195	\$850	Bulk Rate
PSB180NM	0.18	\$103	\$195	\$850	Bulk Rate
PSB300NM	0.3	\$103	\$195	\$850	Bulk Rate
PSB400NM	0.4	\$103	\$195	\$850	Bulk Rate
PSB500NM	0.5	\$103	\$195	\$850	Bulk Rate
PSB800NM	0.8	\$103	\$195	\$850	Bulk Rate
PSB001UM	1	\$120	\$230	\$1,150	Bulk Rate
PSB005UM	5	\$140	\$270	\$1,250	Bulk Rate
PSB007UM	7	\$190	\$375	\$1,600	Bulk Rate
PSB010UM	10	\$240	\$450	\$2,000	Bulk Rate

RED Polystyrene Latex Beads

We also offer colored polystyrene latex beads, which is also in aqueous suspension at 10% solids. You can also contact us if you are looking for a specific color intensity. Product will be shipped within 10% of the diameter you specify.

Catalog No.	Diameter (um)	5 mL	15 mL	100 mL	500 + mL
PSR030NM	0.03	\$103	\$195	\$850	Bulk Rate
PSR050NM	0.05	\$103	\$195	\$850	Bulk Rate
PSR080NM	0.08	\$103	\$195	\$850	Bulk Rate
PSR300NM	0.3	\$103	\$195	\$850	Bulk Rate
PSR400NM	0.4	\$103	\$195	\$850	Bulk Rate
PSR500NM	0.5	\$103	\$195	\$850	Bulk Rate
PSR800NM	0.8	\$103	\$195	\$850	Bulk Rate
PSR001UM	1	\$120	\$230	\$1,150	Bulk Rate
PSR002UM	2	\$120	\$230	\$1,150	Bulk Rate
PSR003UM	3	\$130	\$250	\$1,200	Bulk Rate
PSR005UM	5	\$140	\$270	\$1,250	Bulk Rate
PSR007UM	7	\$190	\$375	\$1,600	Bulk Rate
PSR010UM	10	\$240	\$450	\$2,000	Bulk Rate

Quality manufacturing since 1985

Color Particles

**Blue Carboxylated Latex Beads**-Ideal for Rapid Tests (lateral flow membrane tests)  
We offer colored carboxylated latex beads, which is also in aqueous suspension at 10% solids. You can also contact us if you are looking for a specific color intensity. Product will be shipped within 10% of the diameter you specify. CAB300NM and CAB400NM are most popular among rapid tests manufactures.

Catalog No.	Diameter (um)	5 mL	15 mL	100 mL	500 + mL
CAB150NM	0.15	\$110	\$210	\$850	Bulk Rate
CAB200NM	0.2	\$110	\$210	\$850	Bulk Rate
CAB300NM	0.3	\$110	\$210	\$850	Bulk Rate
CAB400NM	0.4	\$110	\$210	\$850	Bulk Rate
CAB500NM	0.5	\$110	\$210	\$850	Bulk Rate
CAB600NM	0.6	\$110	\$210	\$850	Bulk Rate
CAB700NM	0.7	\$110	\$210	\$850	Bulk Rate
CAB800NM	0.8	\$110	\$210	\$850	Bulk Rate
CAB900NM	0.9	\$110	\$210	\$850	Bulk Rate
CAB001UM	1	\$123	\$235	\$1,150	Bulk Rate
CAB002UM	2	\$123	\$235	\$1,150	Bulk Rate
CAB005UM	5	\$143	\$275	\$1,250	Bulk Rate
CAB007UM	7	\$195	\$380	\$1,600	Bulk Rate
CAB010UN	10	\$240	\$455	\$2,000	Bulk Rate

Red Carboxylated Latex Particles

We also offer colored carboxylated latex beads, which is also in aqueous suspension at 10% solids. You can also contact us if you are looking for a specific color intensity. Product will be shipped within 10% of the diameter you specify. Our CAR400NM is considered the gold standards for rapid test manufacturers.

Catalog No.	Diameter (um)	5 mL	15 mL	100 mL	500 + mL
CAR040NM	0.04	\$110	\$210	\$850	Bulk Rate
CAR200NM	0.2	\$110	\$210	\$850	Bulk Rate
CAR300NM	0.3	\$110	\$210	\$850	Bulk Rate
CAR400NM	0.4	\$110	\$210	\$850	Bulk Rate
CAR500NM	0.5	\$110	\$210	\$850	Bulk Rate
CAR700NM	0.7	\$110	\$210	\$850	Bulk Rate
CAR800NM	0.8	\$110	\$210	\$850	Bulk Rate
CAR001UM	1	\$123	\$235	\$1,150	Bulk Rate
CAR002UM	2	\$123	\$235	\$1,150	Bulk Rate
CAR005UM	5	\$143	\$275	\$1,250	Bulk Rate
CAR007UM	7	\$195	\$380	\$1,600	Bulk Rate
CAR010UM	10	\$240	\$455	\$2,000	Bulk Rate

Surfactant free or near surfactant free particles

Surfactant free or near Surfactant free Carboxylated Polystyrene Particles

This product line is designed for application that surfactant may cause interference. For diameter larger than 0.7 um latex particles are synthesized without surfactant and for particles smaller than 0.7um, there is less than 0.02% surfactant present. These are also water suspension at 10% solid, and the average diameter will be within 10% catalog number specified.

Catalog No.	Diameter (um)	5 mL	15 mL	100 mL	500 + mL
SFCA050NM	0.05	\$90	\$170	\$520	Bulk Rate
SFCA090NM	0.09	\$90	\$170	\$520	Bulk Rate
SFCA150NM	0.15	\$90	\$170	\$520	Bulk Rate
SFCA200NM	0.2	\$90	\$170	\$520	Bulk Rate
SFCA300NM	0.3	\$90	\$170	\$520	Bulk Rate
SFCA400NM	0.4	\$90	\$170	\$520	Bulk Rate
SFCA700NM	0.7	\$90	\$170	\$520	Bulk Rate
SFCA900NM	0.9	\$90	\$170	\$520	Bulk Rate
SFCA002UM	2	\$93	\$170	\$575	Bulk Rate

Surfactant free or near Surfactant free Polystyrene Particles

This product line is designed for application that surfactant may cause interference. For diameter larger than 0.7u,latex , particles are synthesized without surfactant and for particles smaller than 0.7um, there is less than 0.02% surfactant present. These are also water suspension at 10% solid, and the average diameter will be within 10% catalog number specified.

Catalog No.	Diameter (um)	5 mL	15 mL	100 mL	500 + mL
SF060NM	0.06	\$85	\$160	\$520	Bulk Rate
SF080NM	0.08	\$85	\$160	\$520	Bulk Rate
SF090NM	0.09	\$85	\$160	\$520	Bulk Rate
SF120NM	0.12	\$85	\$160	\$520	Bulk Rate
SF150NM	0.15	\$85	\$160	\$520	Bulk Rate
SF300NM	0.3	\$85	\$160	\$520	Bulk Rate
SF400NM	0.4	\$85	\$160	\$520	Bulk Rate
SF500NM	0.5	\$85	\$160	\$520	Bulk Rate
SF600NM	0.6	\$85	\$160	\$520	Bulk Rate
SF700NM	0.7	\$85	\$160	\$520	Bulk Rate
SF800NM	0.8	\$85	\$160	\$520	Bulk Rate
SF900NM	0.9	\$85	\$160	\$520	Bulk Rate
SF001UM	1.5	\$90	\$165	\$575	Bulk Rate

Model Nonspherical Latex Particles

Our model nonspherical Latex particles come in of shapes of peanut or snowman or pear. Those Monosized nonspherical polystyrene particles offer unique physical properties especially in light scattering and fluid dynamic behavior. The entire population per batch is identical in size and shape. These are also in aqueous suspension at 10% solids.and length specified here is the length of the longest portion of the particles.

Catalog No.	Diameter (um)	5 mL	15 mL	100 mL	500 + mL
PNT003UM	3	\$130	\$260	\$1,050	Bulk Rate
PNT004UM	4	\$130	\$260	\$1,050	Bulk Rate
PNT005UM	5	\$170	\$330	\$1,550	Bulk Rate
PNT006UM	6	\$170	\$330	\$1,550	Bulk Rate
PNT007UM	7	\$170	\$330	\$1,550	Bulk Rate
PNT008UM	8	\$185	\$400	\$1,700	Bulk Rate
PNT010UM	10	\$220	\$430	\$1,950	Bulk Rate
PNT012UM	12	\$220	\$430	\$1,950	Bulk Rate
PNT015UM	15	\$220	\$500	\$2,250	Bulk Rate

## Fluorescent Latex Particles

### Fluorescent Polystyrene Latex Beads

Our Fluorescent PS latex beads comes in four (4) different colors; Red, Yellow-green, Blue, and Orange. All of them come in aqueous suspension at 2.5% solids. All our fluorescent particle's size are traceable to NIST size standards. Price: \$65/5 ml, \$130/15 ml, \$400/100 ml, bulk rate available.

Nominal Diameter (um)	Red	Yellow-green	Blue	Orange
	Catalog No.	Catalog No.	Catalog No.	Catalog No.
0.05	PSFR050NM	PSYF050NM	PSFB050NM	PSOF050NM
0.1	PSFR100NM	PSYF100NM	PSFB100NM	PSOF100NM
0.3	PSFR300NM	PSYF300NM	PSFB300NM	PSOF300NM
0.5	PSFR500NM	PSYF500NM	PSFB500NM	PSOF500NM
1	PSFR001UM	PSYF001UM	PSFB001UM	PSOF001UM
2	PSFR002UM	PSYF002UM	PSFB002UM	PSOF002UM
5	PSFR005UM	PSYF005UM	PSFB005UM	PSOF005UM
7	PSFR007UM	PSYF007UM	PSFB007UM	PSOF007UM
10	PSFR010UM	PSYF010UM	PSFB010UM	PSOF010UM

### Fluorescent Carboxylated Latex Beads

Our Fluorescent CA latex beads comes in four (4) different colors; Red, Yellow-green, Blue, and Orange. All of them come in aqueous suspension at 2.5% solids. All our fluorescent particle's size are traceable to NIST size standards. Price: \$70/5 ml, \$135/15 ml, \$400/100 ml, bulk rate available.

Nominal Diameter (um)	Red	Yellow-green	Blue	Orange
	Catalog No.	Catalog No.	Catalog No.	Catalog No.
0.05	CAFR050NM	CAYF050NM	CAFB050NM	CAOF050NM
0.1	CAFR100NM	CAYF100NM	CAFB100NM	CAOF100NM
0.3	CAFR300NM	CAYF300NM	CAFB300NM	CAOF300NM
0.5	CAFR500NM	CAYF500NM	CAFB500NM	CAOF500NM
1	CAFR001UM	CAYF001UM	CAFB001UM	CAOF001UM
2	CAFR002UM	CAYF002UM	CAFB002UM	CAOF002UM
5	CAFR005UM	CAYF005UM	CAFB005UM	CAOF005UM
7	CAFR007UM	CAYF007UM	CAFB007UM	CAOF007UM
10	CAFR010UM	CAYF010UM	CAFB010UM	CAOF010UM

### Fluorescent Aminated Latex Beads

Our Fluorescent AM latex beads comes in four (4) different colors; Red, Yellow-green, Blue, and Orange. All of them come in aqueous suspension at 2.5% solids. All our fluorescent particle's size are traceable to NIST size standards. Price: \$70/5 ml, \$140/15 ml, \$400/100 ml, bulk rate available.

Nominal Diameter (um)	Red	Yellow-green	Blue	Orange
	Catalog No.	Catalog No.	Catalog No.	Catalog No.
0.05	AMFR050NM	AMYP050NM	AMFB050NM	AMOF050NM
0.1	AMFR100NM	AMYP100NM	AMFB100NM	AMOF100NM
0.3	AMFR300NM	AMYP300NM	AMFB300NM	AMOF300NM
0.5	AMFR500NM	AMYP500NM	AMFB500NM	AMOF500NM
1	AMFR001UM	AMYP001UM	AMFB001UM	AMOF001UM
2	AMFR002UM	AMYP002UM	AMFB002UM	AMOF002UM

### How to Order:

Telephone: 1-626-584-9446 Fax: 1-626-795-0797

www.magsphere.com

1993 E. Locust Street, Pasadena, CA 91107

## VBC/PS Latex Beads

VBC/PS Latex Particles are also known as polychloromethylstyrene.

VBC allows for direct covalent bondings. Comes in aqueous suspension at 10% solids, and will be shipped within 10% of the diameter you specify.

Catalog No.	Diameter (um)	5 mL	15 mL	100 mL	500 + mL
VB050NM	0.05	\$90	\$160	\$580	Bulk Rate
VB070NM	0.07	\$90	\$160	\$580	Bulk Rate
VB080NM	0.08	\$90	\$160	\$580	Bulk Rate
VB200NM	0.2	\$90	\$160	\$580	Bulk Rate
VB400NM	0.4	\$90	\$160	\$580	Bulk Rate
VB001UM	1	\$103	\$180	\$650	Bulk Rate
VB003UM	3	\$125	\$250	\$850	Bulk Rate
VB005UM	5	\$180	\$350	\$1,550	Bulk Rate

## Magnetic Particles

### Magnetic Microparticles (Magnetic Beads)

Our magnetic particles, or magnetic beads come in aqueous suspension of 2.5% solids. We have two (2) superparamagnetic, which has Polystyrene based, or Carboxylated PS based. Product will be shipped within 10% of the diameter you specify. We offer in 10 mL, and 100 mL quantity. For bulk quantity of 500+ mL please contact Magsphere Inc.

### Polystyrene Based Magnetic Particle

Catalog No.	Size (Nominal Diameter)	10 mL	100 mL
MPS400NM	average 400nm	\$100	\$750
MPS800NM	average 800nm	\$100	\$750
MPS2UM	average 2um	\$100	\$750
MPS5UM	average 5um	\$100	\$750

### Carboxylated based magnetic Particles

Catalog No.	Size (Nominal Diameter)	10 mL	100 mL
MCA400NM	average 400nm	\$120	\$850
MCA800NM	average 800nm	\$120	\$850
MCA2UM	average 2um	\$120	\$850
MCA5UM	average 5um	\$120	\$850

## Pre-activated Latex Particles

### Polyacrolein Latex

This particle can be used for direct covalent binding via their aldehyde groups. Also comes in aqueous suspension at 10% solids.

Catalog No.	Size (Nominal Diameter)	5 mL	15 mL	100 mL	500 + mL
AL125NM	0.125um	\$120	\$230	\$750	Bulk Rate
AL300NM	0.300um	\$120	\$230	\$750	Bulk Rate

### Latex Particles Containing Epoxy Group

Next, we have Latex Particles containing epoxy group. Also comes in aqueous suspension at 10% solids.

Catalog No.	Size (Nominal Diameter)	5 mL	15 mL	100 mL	500 + mL
EP180NM	0.18um	\$120	\$230	\$750	Bulk Rate

### Uniform Hydrazide/PS Latex

The Uniform Hydrazide/PS Latex, comes in aqueous suspension at 10%.

Catalog No.	Size (Nominal Diameter)	5 mL	15 mL	100 mL	500 + mL
HZ350NM	0.350um	\$125	\$240	\$680	Bulk Rate



*Quality manufacturing since 1985*